

## TRAP Assay Set

Tartrate-Resistant Acid Phosphatase (TRAP) is widely recognized as a marker of osteoclast differentiation and is released into the medium during differentiation in osteoclast cultures<sup>1)2)</sup>. This product is a set of reagents to detect TRAP activity in cultured osteoclasts and can be used to evaluate osteoclast differentiation and drugs for osteoporosis. A TRAP substrate solution containing p-nitrophenyl phosphate is added to the culture supernatant or fixed cells, and the p-nitrophenol produced by the reaction of the TRAP enzyme is detected by absorbance at 405 nm in this simple method.

### **1) Composition of reagent sets**

- ①TRAP Assay Substrate Solution (Cat.TRAP-S1, 15mL, 2 bottles)

Citrate buffer (pH 5.5) containing p-nitrophenyl phosphate and sodium tartrate

- ②TRAP Assay Stop Solution (Cat.TRAP-E1, 15mL, 2 bottles)

Sodium hydroxide solution (0.1 mol/L)

### **2) Example of us**

#### **TRAP activity of culture supernatants**

- (1) Inoculate RAW264.7 cells into each well in culture medium (MEMα containing 10% FBS). Add an inducer of osteoclastic differentiation, such as RANKL (100 ng/mL) and the test substances to be evaluated.  
96well : 2,000cells/0.2mL/well  
48well : 5,000cells/0.5mL/well  
24well : 10,000cells/1mL/well
- (2) After 72 or 96 hours of incubation, collect 50 µL of culture supernatant in a 96-well plate for absorbance measurement.
- (3) Add 50 µL of TRAP Assay Substrate Solution to the supernatant of each well, mix and react at 37°C for 30 min to 1 h (light-shielded).
- (4) Add 50 µL of TRAP Assay Stop Solution to each well and stir (colour develops at this point).
- (5) Measure absorbance at 405 nm.

#### **TRAP activity of cells**

- (1) After cultivation, remove the medium from each well of the plate, add 10% neutral buffered formalin solution and fix the cells for 10 min at room temperature.
- (2) Remove the formalin solution and wash the cells in each well twice with purified water (0.2-1 mL).
- (3) Add the following volumes of TRAP Assay Substrate Solution and react at 37°C for 30 min to 1 h (light-shielded).  
96-well plate: 100 µL/well  
48-well plate: 250 µL/well  
24-well plate: 500 µL/well
- (4) Add the same volume of TRAP Assay Stop Solution as Substrate Solution (colour develops at this point) and stir gently.
- (5) After collecting 100 µL of the reaction solution into a 96-well plate, measure the absorbance at 405 nm.

### **3) Assay precautions**

- (1) When using this reagent, always wear protective equipment such as glasses and gloves and take care to avoid contact with the human body.
- (2) After thawing this reagent, stir it well before use.
- (3) When storing residual reagents after measurement use, they should be frozen (-20°C). It has been confirmed that there are no quality problems after refreezing and thawing three times.
- (4) When reacting with TRAP substrates, shield from light and adjust the reaction time to suit your culture system. Colour develops when Stop Solution is added.
- (5) Reactions can be carried out at room temperature, but the overall absorbance will be lower than at 37 °C.
- (6) Can also be used for culture supernatants containing phenol red.
- (7) When using calcium phosphate coated plates (Bone Resorption Assay Plate) to measure TRAP activity of cells, use the supernatant after centrifugation at 2,000 x g for 3 min, as the addition of Stop Solution causes cloudiness.

4) **References**

1. Alatalo SL, Clin Chem, 2000, 46(11):1751-1754.  
2. Lv Y, Exp Ther Med, 2015, 9(1):143-146.

5) **Expected Results**




Fig.1 RAW264.7 cells (48-well plate, 10% FBS/MEMα) were incubated with RANKL for 96 h upon seeding, and the culture supernatant and cells were collected. A: 50 μL of Substrate Solution to 50 μL of culture supernatant, B: 250 μL of Substrate Solution was added to the fixes cells, and allowed to react for 1 h at 37 °C. After adding Stop Solution, absorbance at 405 nm was measured. (Mean±SD, n=4)




Fig.2 RAW264.7 cells (96-well plate, 10% FBS/MEMα) were incubated with RANKL(100ng/mL) and various concentrations of Alendronate for 72 h upon seeding, and the culture supernatant and cells were collected. A: 50 μL of Substrate Solution to 50 μL of culture supernatant, B: 100 μL of Substrate Solution was added to the fixed cells, and allowed to react for 1 h at 37 °C. After adding Stop Solution, absorbance at 405 nm was measured. (Mean±SD, n=4)

| Catalog Number | Product Name   | Specification                       | Price   | Storage |
|----------------|----------------|-------------------------------------|---------|---------|
| TRAP-SET       | TRAP Assay Set | 1 set<br>(TRAP-S1, TRAP-E1: 2 each) | ¥21,000 | -20℃    |